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The image shows two overlapping screenshots from the Mascot web interface. The background screenshot is the 'Mascot Utilities' page, which has a sidebar with links: 'Database Status', 'Search Log', 'Configuration Editor' (highlighted with a red box), and 'Database Manager'. The foreground screenshot is the 'Mascot Configuration' page, which displays a table of configuration options. The 'Database Manager' link at the bottom of this table is also highlighted with a red box.

Amino Acids	Amino Acid Data
Modifications	Modification definitions
Symbols	Symbols used in chemical formulae
Linkers	Linker definitions
Enzymes	Enzyme definitions
Instruments	Fragmentation Rules
Quantitation	Quantitation Methods
Crosslinking	Crosslinking Methods
Configuration Options	Global Options in mascot.dat
Database Manager	Sequence databases, Parse Rules

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```

Name      = PRIDE_Contaminants  Family   = C:/inetpub/mascot/sequence/PRIDE_Contaminants/current/PRIDE_Contaminants
Filename  = PRIDE_Contaminants_20160908.msp Pathname = C:/inetpub/mascot/sequence/PRIDE_Contaminants/current/PRIDE_Contaminants_20160908.msp
Status    = In use
State Time = Tue Apr 13 09:33:38 # searches = 0
Mem mapped = NO Request to mem map = NO Request unmap = NO Mem locked = NO
Number of threads = 1 Current = YES Type = Spectral library

Name      = SwissProt          Family   = C:/inetpub/mascot/sequence/SwissProt/current/SwissProt_*.fasta
Filename  = SwissProt_2020_05.fasta Pathname = C:/inetpub/mascot/sequence/SwissProt/current/SwissProt_2020_05.fasta
Status    = In use
State Time = Tue Apr 13 09:33:38 # searches = 0
Mem mapped = YES Request to mem map = YES Request unmap = NO Mem locked = NO
Number of threads = -1 Current = YES Type = Amino acid

Name      = UP5640_H_sapiens   Family   = C:/inetpub/mascot/sequence/UP5640_H_sapiens/current/UP5640_H_sapiens
Filename  = UP5640_H_sapiens_20201007.fasta Pathname = C:/inetpub/mascot/sequence/UP5640_H_sapiens/current/UP5640_H_sapiens_20201007.fasta
Status    = In use
State Time = Tue Apr 13 09:33:38 # searches = 0
Mem mapped = YES Request to mem map = YES Request unmap = NO Mem locked = NO
Number of threads = -1 Current = YES Type = Amino acid

Name      = OriginalDB        Family   = C:/inetpub/mascot/sequence/OriginalDB/current/OriginalDB_*.fasta
Filename  = OriginalDB_20210412.fasta Pathname = C:/inetpub/mascot/sequence/OriginalDB/current/OriginalDB_20210412.fasta
Status    = In use
State Time = Tue Apr 13 11:20:04 # searches = 0
Mem mapped = YES Request to mem map = YES Request unmap = NO Mem locked = NO
Number of threads = -1 Current = YES Type = Amino acid

```

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```

Name      = OriginalDB        Family   = C:/INETPUB/MASCOT/sequence/OriginalDB/current/OriginalDB_*.fasta
Filename  = OriginalDB_20210412.fasta Pathname = C:/INETPUB/MASCOT/sequence/OriginalDB/current/OriginalDB_20210412.fasta
Status    = Halted: First test failed. check error log fix fault and retry
State Time = Tue Apr 13 11:20:04 # searches = 0
Mem mapped = NO Request to mem map = YES Request unmap = NO Mem locked = NO
Number of threads = 4 Current = NO

```

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