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
**Mascot Database Manager Configuration Import**

**Note:**  
Please disable other programs and scripts that configure or update sequence databases before continuing.

Database Manager can keep your sequence database configuration up to date and schedule the download of new files. If you choose to use Database Manager, you should not edit `mascot.dat` in future, because your changes will be lost. If you do not wish to use Database Manager now, please [return to the Configuration Editor](#).

Database Manager has two kinds of database definitions:

1. *Predefined database definition*: Automatic configuration by downloading up-to-date information from the Matrix Science website.
2. *Custom database definition*: Manually maintained configuration.

You can choose to keep the existing definition (Keep as Custom) or synchronise with a predefined definition (Synchronise). Synchronisation means the predefined database definition will be applied, and the database files (FASTA, reference, taxonomy) will be updated. If you do not want to update the database files in the future by choosing Synchronise Custom, please choose Keep as Custom.

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Import database configuration by choosing Import.

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
The screenshot shows the Mascot Database Manager interface. The left sidebar has a menu with 'Databases (7)' highlighted. The main content area is titled 'Databases and spectral libraries' and contains a table of database entries. Below the table, there is a note about the latest predefined definitions files and a 'Refresh' button.

Name	Mode	Type	Status	Latest task
contaminants	predefined	AA	In use	Update succeeded (view log)
MyLibrary01	custom	SL	In use	Update succeeded (view log)
neXtProt	predefined	AA	In use	Update succeeded (view log)
NIST_Human_HCD	predefined	SL	In use	Update succeeded (view log)
PRIDE_Contaminants	predefined	SL	In use	Update succeeded (view log)
SwissProt	predefined	AA	In use	Update succeeded (view log)
filters2	custom	SL	Offline	Configuration error

Latest predefined definitions files are from Thu Jul 6 14:23:09 2017 (databases\_1.xml) and Mon May 22 14:16:08 2017 (libraries\_1.xml).  
Full database status is available on [the database status page](#).

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Name	Mode	Type	Status			Latest task
<a href="#">contaminants</a>	predefined	AA	In use	Get new files	Deactivate	Update succeeded ( <a href="#">view log</a> )
<a href="#">MyLibrary01</a>	custom	SL	In use		Deactivate	Update succeeded ( <a href="#">view log</a> )
<a href="#">neXtProt</a>	predefined	AA	In use	Get new files	Deactivate	Update succeeded ( <a href="#">view log</a> )
<a href="#">NIST_Human_HCD</a>	predefined	SL	In use	Get new files	Deactivate	Update succeeded ( <a href="#">view log</a> )
<a href="#">PRIDE_Contaminants</a>	predefined	SL	In use	Get new files	Deactivate	
<a href="#">SwissProt</a>	predefined	AA	In use	Get new files	Deactivate	
<a href="#">rs2</a>	custom	SL	Offline		Activate	

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Matrix Science - Home x Mascot configuration x Mascot Database Ma... x

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**Most recent finished task**  
[Wed May 9 03:17:26 2012] (success) No new FASTA file available; database SwissProt is up to date.

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
The image consists of three overlapping screenshots illustrating the steps to edit the update schedule for a database.

**Step 1:** The top screenshot shows the 'Databases' management page. A table lists databases with columns for Name, Status, New files, and Tasks. The 'SwissProt' database is highlighted with a red box, and a yellow callout box with the text '① DB名をクリック' (Click the DB name) points to it.

**Step 2:** The middle screenshot shows the 'Scheduled updates' section for the SwissProt database. The 'Edit schedule' button is highlighted with a red box, and a yellow callout box with the text '② 「Edit schedule」' (Click 'Edit schedule') points to it.

**Step 3:** The bottom screenshot shows the 'Database configuration: SwissProt' page. Under the 'Scheduled updates' section, the 'Monthly' option is selected with a radio button. The 'Monthly' configuration is set to 'on the 25th at 00:45', which is highlighted with a red box. The 'Save' button is also highlighted with a red box, and a yellow callout box with the text '③ 定期更新のスケジュール設定後「Save」' (After setting the update schedule, click 'Save') points to it.

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
Title	Started	Progress
Download files for Human_EST dbman_download.pl Human_EST	Tue Feb 28 15:04:01 2017	[0.0%]

Process end time	Command	Last message
[Mon Feb 27 18:26:09 2017]	dbman_download.pl PRIDE_Contaminants	(success) 'PRIDE_Contaminants' is up to date.
[Mon Feb 27 18:25:33 2017]	dbman_download.pl SwissProt	(success) 'SwissProt' successfully updated.
[Mon Feb 27 18:18:46 2017]	Predefined definitions file download	(success) Latest predefined definitions file saved to disk.

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The screenshot shows the Mascot Database Manager web interface. The left sidebar contains a menu with the following items: Database Manager, Databases (7), Parse rules (11), Scheduled updates (0), Running tasks (0), **Settings** (highlighted with a red box), Fasta, Enable predefined definition, Synchronise custom definitions, Create new, Library, Enable predefined definition, Synchronise custom definitions, Create new, Spectral library filters.

The main content area is titled "Databases and spectral libraries" and contains a table with the following data:

Name	Mode	Type	Status	Latest task
contaminants	predefined	AA	In use	Update succeeded (view log)
MyLibrary01	custom	SL	In use	Update succeeded (view log)
neXtProt	predefined	AA	In use	Update succeeded (view log)
NIST_Human_HCD	predefined	SL	In use	Update succeeded (view log)
PRIDE_Contaminants	predefined	SL	In use	Update succeeded (view log)
SwissProt	predefined	AA	In use	Update succeeded (view log)
filters2	custom	SL	Offline	Configuration error

Below the table, there is a note: "Latest predefined definitions files are from Thu Jul 6 14:23:09 2017 (databases\_1.xml) and Mon May 22 14:16:08 2017 (libranes\_1.xml). Full database status is available on the database status page." A "Refresh" button is located at the bottom of the table area.

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**Internet access settings**

**Allow Internet access?**

Database Manager may access the Internet for two purposes:

1. Download the latest configuration information from the Matrix Science website.
2. Download new database files.

If the Mascot Server machine has no Internet connection, or if you do not wish Database Manager to access the Internet, this behaviour can be disabled.

Internet access is **enabled**.

Changing this setting will take effect immediately. Update tasks that are currently running will be allowed to continue until they are done.

**Allow external full-text reports?**

By default, if a database definition has an external full-text report source (an external URL), Mascot report scripts will download and display the extra data when protein details are viewed. This behaviour can be disabled.

External full-text report access is **enabled**.

Changing this setting will take effect immediately for all active and inactive database definitions.

**Proxy settings**

You may need to define proxy settings to enable Database Manager to access the Internet, depending on your local network and Internet connection.

Proxy settings defined in `mascot.dat`:

- HTTP proxy strategy: *Auto*
- FTP proxy strategy: *direct connection*

URL test **succeeds**; Database Manager can access [http://www.matrixscience.com/downloads/databases\\_1.xml](http://www.matrixscience.com/downloads/databases_1.xml) from the Mascot server machine. (FTP proxy settings not tested.)

**Sequence directory settings**

Sequence databases are usually stored in a single sequence directory. The default path is `C:/inetpub/mascot/sequence`, but this can be changed below.

Note that you can always change the sequence directory for individual database definitions when creating a new database definition.

Sequence directory:

(If you leave the value empty, it will be reset to the default path.)

**Spectral library settings**

All spectral libraries must have a reference database (FASTA). The reference database is chosen when the library is created (custom definition) or enabled (predefined definition).

If the reference database does not exist or is not a protein sequence database (AA), then Database Manager will offer SwissProt as the default choice.

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**Allow external full-text reports?**

By default, if a database definition has an external full-text report source (an external URL), Mascot report scripts will download and display the extra data when processing. This behaviour can be disabled.

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

(If you leave this field blank, the path will be reset to the default path.)

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The screenshot shows the Mascot Database Manager interface. The left sidebar has a 'Fasta' menu item highlighted with a red box. The main content area is titled 'Databases and spectral libraries' and contains a table with the following data:

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Below the table, there is a 'Refresh' button and a message: 'Latest predefined definitions files are from Thu Jul 6 14:23:09 2017 (databases\_1.xml) and Mon May 22 14:16:08 2017 (libraries\_1.xml). Full database status is available on the database status page.'

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
The screenshot shows the 'Mascot Database Manager' interface. The main heading is 'Enable predefined database definition'. Below the heading, there is a paragraph explaining that predefined definitions are configuration entries for commonly used databases. A table lists various databases with 'Enable' buttons. Two yellow callout boxes provide instructions: ① 'Enable predefined definition' をクリックし and ② 使用したいデータベースの 'Enable' ボタンを押す. The 'Enable' button for 'neXtProt' is highlighted with a red box.


Name	Enable	Notes
contaminants	Enable	Could be synchronised with NCBINr.
cRAP	Enable	Already set up as cRAP.
Environmental_EST	Enable	
Fungi_EST	Enable	
Human_EST	Enable	Already set up as Human_EST.
Invertebrates_EST	Enable	
IPI_arabidopsis	Enable	
IPI_bovine	Enable	
IPI_chicken	Enable	
IPI_human	Enable	Already set up as IPI_human.
IPI_mouse	Enable	
IPI_rat	Enable	
IPI_zebrafish	Enable	
Mammals_EST	Enable	
Mus_EST	Enable	
NCBIprot	Enable	Already set up as NCBIprot.
neXtProt	Enable	
Plants_EST	Enable	

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### Database Manager tasks

Database Manager processes scheduled tasks one at a time in the background as long as Mascot (ms-monitor.exe) is running.

Task processing is currently **enabled**. Task queue and currently running tasks are shown below. ▶ **Task process control**

#	Queued task	Added	
1	Download files for IPI_bovine	Wed May 30 10:57:55	<input type="button" value="Cancel"/>

Title	Started	Progress	
Download files for Environmental_EST dbman_download.pl Environmental_EST	Wed May 30 10:33:01 2012	[12.0%] Downloading <ftp://ftp.ebi.ac.uk/pub/databases/embl/misc/acc_to_taxid.mapping.txt.gz> (try 1 out of 5): 210 kB/s (296 MB/568 MB, 22m 5s remaining)	<input type="button" value="Cancel"/>

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```

Name      = IPI_human          Family   = C:/inetpub/mascot/sequence/IPI_human/current/IPI_human_*.fa
Filename  = IPI_human_3.87x.fasta Pathname = C:/inetpub/mascot/sequence/IPI_human/current/IPI_human_3.8
Status    = In use           Statistics
State Time = Tue May 29 10:59:57 # searches = 0
Mem mapped = YES Request to mem map = YES Request unmap = NO Mem locked = NO
Number of threads = 4 Current = YES

```

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```

Name      = IPI_mouse          Family   = C:/inetpub/mascot/sequence/IPI_mouse/current/IPI_mouse_*.fa
Filename  = IPI_mouse_3.87.fasta Pathname = C:/inetpub/mascot/sequence/IPI_mouse/current/IPI_mouse_3.87
Status    = In use           Statistics
State Time = Tue May 29 10:59:57 # searches = 0
Mem mapped = YES Request to mem map = YES Request unmap = NO Mem locked = NO
Number of threads = -1 Current = YES

```

```

Name      = NCBIInr          Family   = C:/inetpub/mascot/sequence/NCBIInr/current/NCBIInr_*.fasta

```

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The screenshot shows a web browser window with the URL `t5600-1/mascot/x-cgi/db_manager.pl?sub=dbs.merge;dbs.new.def-class=db`. The page title is "Synchronise database definitions". The left sidebar has a menu with "Synchronise custom definitions" highlighted in red. The main content area has the following text:

The following custom definitions match predefined definitions and could be synchronised.

The advantage of synchronisation is automatic updates. As long as the Mascot server machine is connected to the Internet, configuration of the synchronised entry will be automatically kept up to date with the Matrix Science website, and the database files (FASTA, reference, taxonomy) can be automatically downloaded at the frequency of your choice.

Synchronisation will keep local settings fixed, such as directory names. Database files and parse rules will be updated to the latest version. The database can stay online while synchronisation is in progress, but its configuration cannot be changed.

	Action	Match quality
NCBIInr	<input checked="" type="radio"/> Keep <input type="radio"/> Synchronise with contaminants	poor ▶ Differences

At the bottom left of the main content area, there is a "Do it" button.

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The screenshot shows the Mascot Database Manager interface. On the left, the 'Database Manager' sidebar has a 'Create new' button highlighted with a red box. The main content area is titled 'New database definition'. The 'Database name' field contains 'Uniprot\_Human'. Below this, there are radio buttons for 'Custom' and 'Use predefined definition template'. The 'Use predefined definition template' option is selected and highlighted with a red box, with a yellow callout box pointing to it containing the text '「Use predefined definition template」'. Below this, a dropdown menu shows 'UniProt\_proteome\_template'. At the bottom, a 'Next' button is highlighted with a red box.

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
Main file URL  
Delete original file after updating  
Version file URL or path to source file on Mascot Server hard disk (?)  
Delete original file after updating  
Reference file URL or path to source file on Mascot Server hard disk (?)  
Delete original file after updating  
The original file can only be deleted if it resides on the Mascot Server hard disk and Database Manager has sufficient permissions in the source directory.  
Cancel Save

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The screenshot shows the Mascot Database Manager interface for the 'Uniprot\_Human' database. The 'FASTA file to upload' section is highlighted, showing a file selection button circled in red. Below it, the 'Upload' button is also circled in red. Two yellow callout boxes provide instructions: '① FASTA ファイルを選び、' and '② Upload ボタンを押す'.

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**Filename pattern**: Uniprot\_Human\_\*.fasta

Please,

1. copy the FASTA file to C:/inetpub/mascot/sequence/Uniprot\_Human/current ;
2. rename the file to match the filename pattern, e.g. Uniprot\_Human\_20170808.fasta (in the pattern, \* matches any string of characters);
3. and refresh this page to continue.

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The image shows two screenshots from the Mascot database management interface. The top screenshot displays the 'Database: Uniprot\_Human' configuration page. The 'Database status' section indicates the database is 'Offline', and the 'Activate' button is highlighted with a red circle. A yellow callout box contains the text: **ファイル名が認識されたら「Activate」**. The bottom screenshot shows the 'Database Status' page, which lists the status of various databases. The 'Uniprot\_Human' entry is circled in red, and a yellow callout box contains the text: **Database Status 画面で データベース構築確認**.

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The screenshot shows the NCBI Taxonomy browser interface. The search bar contains 'Rattus norvegicus'. The search results show the following information:

**Rattus norvegicus**

Taxonomy ID: 10116  
 Genbank common name: Norway rat  
 Inherited blast name: rodents  
 Rank: species  
 Genetic code: Translation table 1 (Standard)  
 Mitochondrial genetic code: Translation table 2 (Vertebrate Mitochondrial)  
 Other names:  
 common name: brown rat

**Entrez records**

Database name	Subtree links	Direct links
Nucleotide	<a href="#">270,518</a>	<a href="#">270,517</a>
Nucleotide EST	<a href="#">1,103,577</a>	<a href="#">1,103,577</a>
Nucleotide GSS	<a href="#">867,131</a>	<a href="#">867,131</a>
Protein	<a href="#">110,744</a>	<a href="#">110,744</a>
Structure	<a href="#">1,736</a>	<a href="#">1,736</a>
Genome	<a href="#">44</a>	<a href="#">44</a>

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Download selected (0)  
 Download all (68511)  
 Format: FASTA (canonical)  
 Preview first 10  
 Go

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

Name	Mode	Type	Status	Latest task
contaminants	predefined	AA	In use	Update succeeded (view log)
MyLibrary01	custom	SL	In use	Update succeeded (view log)
neXtProt	predefined	AA	In use	Update succeeded (view log)
NIST_Human_HCD	predefined	SL	In use	Update succeeded (view log)
PRIDE_Contaminants	predefined	SL	In use	Update succeeded (view log)
SwissProt	predefined	AA	In use	Update succeeded (view log)
filters2	custom	SL	Offline	Configuration error
			Offline	No local FASTA file

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**Database direct**  
C:/inetpub/mascot

**Filename patter**  
OriginalDB2\_\*.fas

Database files must be activated. You have the following options:

- Download from remote URL or copy from Mascot server hard disk
- Upload file using web browser
- Copy file manually

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Name  
OriginalDB2

Database type  
Amino acid (protein database)

Database directory  
C:/inetpub/mascot/sequence/OriginalDB2/current

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The screenshot shows the Mascot Database Manager interface. The main content area displays the configuration for 'Database: OriginalDB2'. Below this, the 'Upload FASTA file' section is visible. The 'FASTA file to upload' section has a '参照...' button circled in red. A yellow callout box with red text says '① FASTA ファイルを選び、'. Below the 'FASTA file to upload' section, the 'Upload' button is circled in red. Another yellow callout box with red text says '② Upload ボタンを押す'.

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The screenshot shows the Mascot web interface for a database named 'OriginalDB2'. The interface is divided into several sections:

- Database Manager:** Includes links for Databases (10), Parse rules (15), Scheduled updates (0), Running tasks (0), and Settings.
- Fasta:** Includes links for Enable predefined definition, Synchronise custom definitions, and Create new.
- Library:** Includes links for Enable predefined definition, Synchronise custom definitions, Create new, and Spectral library filters.

The main configuration area for 'Database: OriginalDB2' includes:

- Name:** OriginalDB2
- Database type:** Amino acid (protein database)
- Database directory:** C:/inetpub/mascot/sequence/OriginalDB2/current (highlighted with a red box and a yellow callout: ① 以下で指定されたフォルダに)
- Filename pattern:** OriginalDB2\_\*.fasta

The 'Copy files manually' section includes:

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
An inset window shows a file explorer view of the 'current' directory, containing a file named 'OriginalDB2\_2\_20170810.fasta'.

At the bottom of the page, there are 'Previous' and 'Refresh' buttons.

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The screenshot shows a web browser window with the URL `localhost/mascot/x-cgi/`. The page title is "Database: OriginalDB2". On the left, there is a navigation menu with sections: "Database Manager" (Databases (10), Parse rules (15), Scheduled updates (0), Running tasks (0), Settings), "Fasta" (Enable predefined definition, Synchronise custom definitions, Create new), and "Library" (Enable predefined definition, Synchronise custom definitions, Create new, Spectral library filters). The main content area displays the following information:

- Copy** **Delete**
- Name**  
*OriginalDB2*
- Database type**  
Amino acid (protein database)
- Database directory**  
C:/inetpub/mascot/sequence/OriginalDB2/current
- Filename pattern**  
OriginalDB2\_\*.fasta
- Files matching OriginalDB2\_\*.fasta**  
OriginalDB2\_20170810.fasta (361 B)
- Upload files**
- Database status**  
Offline  
**Activate** (button circled in red)
- Scheduled updates**  
(no schedules defined)  
**Update now**
- Main database file URL has not been defined; database updates are disabled.
- Show configuration details**

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The screenshot shows the Mascot Database Status page with the following information:

```
State Time = Wed Aug 09 19:04:14 # searches = 0
Mem mapped = NO Request to mem map = NO Request unmap = NO Mem locked = NO
Number of threads = 1 Current = YES Type = Spectral library
```

Database 1:

```
Name = PRIDE_Contaminants Family = C:/inetpub/mascot/sequence/PRIDE_Contaminants/current/PRIDE_C
Filename = PRIDE_Contaminants_20160908.msp Pathname = C:/inetpub/mascot/sequence/PRIDE_Contaminants/curr
Status = In use Statistics Recompress file
State Time = Wed Aug 09 19:04:14 # searches = 0
Mem mapped = NO Request to mem map = NO Request unmap = NO Mem locked = NO
Number of threads = 1 Current = YES Type = Spectral library
```

Database 2:

```
Name = SwissProt Family = C:/inetpub/mascot/sequence/SwissProt/current/SwissProt *.fasta
Filename = SwissProt_2017_05.fasta Pathname = C:/inetpub/mascot/sequence/SwissProt/current/SwissProt_201
Status = In use Statistics Unidentified taxonomy Recompress file
State Time = Wed Aug 09 19:04:14 # searches = 0
Mem mapped = YES Request to mem map = YES Request unmap = NO Mem locked = NO
Number of threads = -1 Current = YES Type = Amino acid
```

Database 3:

```
Name = OriginalDB Family = C:/inetpub/mascot/sequence/OriginalDB/current/OriginalDB_*.fa
Filename = OriginalDB_20170810.fasta Pathname = C:/inetpub/mascot/sequence/OriginalDB/current/OriginalDB
Status = In use Statistics Recompress file
State Time = Mon Aug 14 11:43:22 # searches = 0
Mem mapped = YES Request to mem map = YES Request unmap = NO Mem locked = NO
Number of threads = -1 Current = YES Type = Amino acid
```

Database 4 (highlighted):

```
Name = OriginalDB2 Family = C:/inetpub/mascot/sequence/OriginalDB2/current/OriginalDB2_*.
Filename = OriginalDB2_20170810.fasta Pathname = C:/inetpub/mascot/sequence/OriginalDB2/current/Original
Status = In use Statistics Recompress file
State Time = Mon Aug 14 14:04:24 # searches = 0
Mem mapped = YES Request to mem map = YES Request unmap = NO Mem locked = NO
Number of threads = -1 Current = YES Type = Amino acid
```

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

The screenshot shows the 'Custom definition' page in the Mascot Database Manager. The page is titled 'Custom definition' and has a sidebar on the left with navigation links: Database Manager, Databases (8), Parse rules (13), Scheduled updates (0), Running tasks (0), Settings, Fasta, Library, and Spectral library filters. The main content area is for configuring a new database. The 'Database name' is 'OriginalDB'. The 'Custom definition type' is 'Amino acid (protein database)'. The 'Base directory' is 'C:/inetpub/mascot/sequence'. There are three radio button options for file acquisition: 'Download or copy files automatically', 'Upload or copy files manually', and 'Reference file URL or path to source file on Mascot Server hard disk (optional)'. The 'Upload or copy files manually' option is selected. The 'Create' button is highlighted with a red box. Several yellow callout boxes provide Japanese annotations: '「Sequence directory」配列データベースファイルの置き場所' points to the base directory field; '「Download database files」自動更新プログラムにおけるファイルの取得先' points to the 'Download or copy files automatically' option; '「Upload or Copy files manually」上記設定にて自動ダウンロードしない場合に選択' points to the selected 'Upload or copy files manually' option; and '入力、選択後「Create」' points to the 'Create' button.

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**Database directory**  
C:/inetpub/mascot/sequence/OriginalDB/current

**Filename pattern**  
OriginalDB\_\*.fasta

Database files must be present before database configuration can continue. You have the following options.

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Name  
OriginalDB

Database type  
Amino acid (protein database)

Database directory  
C:/inetpub/mascot/sequence/OriginalDB/current

Filename pattern  
OriginalDB\_\*.fasta

### Download files from URL or copy from Mascot Server disk

The FASTA file can also be downloaded from a remote URL or copied from the Mascot Server hard disk.

Set up download URL

Filename pattern  
OriginalDB\_\*.fasta

Edit filename pattern

The downloaded files will be renamed to match the filename pattern.

Previous Start downloading

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### Database configuration: OriginalDB

Main file URL or path to source file on Mascot Server hard disk

Delete original file after updating

Version file URL or path to source file on Mascot Server hard disk (?)

Delete original file after updating

Reference file URL or path to source file on Mascot Server hard disk (?)

Delete original file after updating

The original file can only be deleted if it resides on the Mascot Server hard disk and Database Manager has sufficient permissions in the source directory.

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The screenshot shows the Mascot Database Manager interface. The browser address bar is `localhost/mascot/x-cgi/db_manager.pl?dbs.file-source=upload&st`. The main content area is titled "Database: OriginalDB" and includes a sidebar with navigation options like "Database Manager", "Fasta", and "Library". The "FASTA file to upload" section is highlighted with a yellow box and a red circle around the "参照..." button. A yellow callout box with the text "① FASTA ファイルを選び、" points to the "参照..." button. Another yellow callout box with the text "② Upload ボタンを押す" points to the "Upload" button.

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

Accession parse rule (?)  
(none)  
**Choose**

Description parse rule (?)  
(none)  
**Choose**

**Accession parse rule for OriginalDB**

**Raw test data (4 samples)**

- >ID0001
- >ID0002 Description 2
- >ID0003 Shortpeptide
- >ID0004 TEST004

**Suitable parse rules**

Parse rule	Match	Extracted data
<input type="radio"/> >¥(.)*¥	4/4	1. ID0001 ▶ 3 more matches (Warning: contains spaces)
<input type="radio"/> >¥([^\ ])*¥	4/4	1. ID0001 ▶ 3 more matches
<input checked="" type="radio"/> >¥([^\ ])*¥	4/4	1. ID0001 ▶ 3 more matches
<input type="radio"/> >¥([^\ ])*¥	4/4	1. ID0001 ▶ 3 more matches
<input type="radio"/> ¥(.)*¥	4/4	▶ 3 more matches
<input type="radio"/> ¥([^\ ])*¥	4/4	▶ 3 more matches

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**Description parse rule (?)**  
(none)  
Choose  
Cancel Next  
Accession and description parse next step.

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
**Database configuration: OriginalDB (step 2/2)**

**Taxonomy source**

- None (?)
- FASTA file (?)  
All human with TaxID 9606

**Sequence report source**

FASTA file

**Full-text report source**

- None (?)
- External source (?)

**Full-text report source**

- None (?)
- External source (?)
  - Current setting  
URL template  
  
Parse rule  
(none)
  - Copy from *Environmental\_EST* (and similar entries)  
<http://www.ebi.ac.uk/cgi-bin/embifetch?id=#ACCESSION#>  
(example)
  - Copy from *Trembl\_LAC* (and similar entries)  
<http://www.uniprot.org/uniprot/#ACCESSION#.txt>  
(example)
  - Copy from *EST\_human* (or *EST\_mouse* or *EST\_others*)  
<https://eutils.ncbi.nlm.nih.gov/entrez/efetch.fcgi?rettype=gb&retmode=text&db=nucleotide&tool=mascot&email=support@matrixscience.com&id=#ACCESSION#>  
(example)
  - Copy from *NCBIInr* (or *NCBIprot*)  
<https://eutils.ncbi.nlm.nih.gov/entrez/efetch.fcgi?rettype=gb&retmode=text&db=protein&tool=mascot&email=support@matrixscience.com&id=#ACCESSION#>  
(example)
- Local program (?)  
Command line template  
  
Parse rule  
(none)  
**To choose or edit local program as a full-text report source, please bring the database online.**

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The image shows two screenshots from the Mascot software interface. The top screenshot is the 'Database: OriginalDB' configuration page. It shows various settings like Name, Database type, Directory, and Filename pattern. The 'Database status' is currently 'Offline'. A red box highlights the 'Activate' button, and a yellow callout box with a circled '9' says '「Activate」 Database 構築開始'.

The bottom screenshot is the 'Mascot search status page' showing a list of databases. The 'OriginalDB' entry at the bottom is highlighted with a red box, and its status is 'In use'. A yellow callout box with a circled '10' says '目的データベースの status 項目が「In use」となれば 使用可能'.

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**Database status**  
Offline

Activate

**Scheduled updates**  
(no schedules defined)

Update now

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Mascot search engine | Mascot search status page | Mascot configuration | Mascot Database Manager

localhost/mascot/x-cgi/ms-status.exe

Status = In use [Statistics](#) [Recompress file](#)

State Time = Wed Aug 09 19:04:14 # searches = 0

Mem mapped = NO Request to mem map = NO Request unmap = NO Mem locked = NO

Number of threads = 1 Current = YES Type = Spectral library

Name = PRIDE\_Contaminants\_ Family = C:/inetpub/mascot/sequence/PRIDE\_Contaminants/current/PRIDE\_Contaminants\_\*.fasta

Filename = PRIDE\_Contaminants\_20160908.msp Pathname = C:/inetpub/mascot/sequence/PRIDE\_Contaminants/current/PRIDE\_Contaminants\_20160908.msp

Status = In use [Statistics](#) [Recompress file](#)

State Time = Wed Aug 09 19:04:14 # searches = 0

Mem mapped = NO Request to mem map = NO Request unmap = NO Mem locked = NO

Number of threads = 1 Current = YES Type = Spectral library

Name = SwissProt\_ Family = C:/inetpub/mascot/sequence/SwissProt/current/SwissProt\_2017\_06.fasta

Filename = SwissProt\_2017\_06.fasta Pathname = C:/inetpub/mascot/sequence/SwissProt/current/SwissProt\_2017\_06.fasta

Status = In use [Statistics](#) [Unidentified taxonomy](#) [Recompress file](#)

State Time = Wed Aug 09 19:04:14 # searches = 0

Mem mapped = YES Request to mem map = YES Request unmap = NO Mem locked = NO

Number of threads = -1 Current = YES Type = Amino acid

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Name = OriginalDB\_ Family = C:/inetpub/mascot/sequence/OriginalDB/current/OriginalDB\_\*.fasta

Filename = OriginalDB\_20170810.fasta Pathname = C:/inetpub/mascot/sequence/OriginalDB/current/OriginalDB\_20170810.fasta

Status = In use [Statistics](#) [Recompress file](#)

State Time = Thu Aug 10 11:23:17 # searches = 0

Mem mapped = YES Request to mem map = YES Request unmap = NO Mem locked = NO

Number of threads = -1 Current = YES Type = Amino acid

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The screenshot shows the Mascot Database Manager interface. The main content area is titled 'Databases and spectral libraries' and contains a table of database entries. The left sidebar has a 'Library' menu item highlighted with a red circle. Below the table, there is a note about predefined definitions and a 'Refresh' button.

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

The screenshot shows the Mascot Database Manager interface. The main heading is "Enable predefined library definition". Below the heading, there is a table of predefined libraries. The table has two columns: "Name" and "Enable". The "Enable" column contains buttons for each library. The "Enable predefined definition" link in the left sidebar is highlighted with a red box. The "Enable" button for "NIST\_Mouse\_HCD" is also highlighted with a red box. Two yellow callout boxes provide instructions: ① 「Enable predefined definition」をクリックし (Click 'Enable predefined definition') and ② 使用したいデータベースの「Enable」ボタンを押す (Press the 'Enable' button for the database you want to use).

Name	Enable
NIST_BSA_IonTrap	Enable
NIST_C.elegans_IonTrap	Enable
NIST_Chicken_IonTrap	Enable
NIST_D.rerio_IonTrap	Enable
NIST_Drosophila_IonTrap	Enable
NIST_E.coli_IonTrap	Enable
NIST_HSA_IonTrap	Enable
NIST_Human_HCD	Enable <small>Already set up as NIST_Human_HCD.</small>
NIST_Mouse_HCD	Enable
NIST_Mouse_HCD_iTRAQ	Enable

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The screenshot shows a web browser window with the URL `localhost/mascot/x-cgi/db_manager.php`. The main content area is titled "Database: NIST\_Mouse\_HCD" and includes several sections:

- Database Manager**: A sidebar menu with links for Databases (11), Parse rules (15), Scheduled updates (0), Running tasks (1), and Settings.
- Fasta**: A sidebar menu with links for Enable predefined definition, Synchronise custom definitions, and Create new.
- Library**: A sidebar menu with links for Enable predefined definition, Synchronise custom definitions, Create new, and Spectral library filters.
- Database configuration details**:
  - Name**: `NIST_Mouse_HCD`
  - Using predefined definition**: `NIST_Mouse_HCD`, Configuration last updated Mon May 22 14:16:08 2017.
  - Database type**: Spectral library
  - Database directory**: `C:/inetpub/mascot/sequence/NIST_Mouse_HCD/current`
  - Filename pattern**: `NIST_Mouse_HCD_*.msp`
- About to download or copy files**:
  - Database files need to be downloaded or copied before database configuration can be completed.
  - Download task is in queue position **1**; waiting for downloading to start. See [the task queue](#) for more details. (Download task added Mon Aug 14 16:21:42 2017.)

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The screenshot shows the Mascot Database Manager interface for a database named "NIST\_Human\_HCD\_20170816". The interface includes a sidebar with navigation options like "Database Manager", "Fasta", and "Library". The main content area shows the database name, type (Spectral library), and directory. Below this, there is a section for "Download files from URL or copy from Mascot Server disk" which includes a "FASTA file URL" field with a value starting with "ftp://chemdata.nist.gov/download/pep...", a red-bordered "Edit download URLs" button, and a "Filename pattern" field with a value "NIST\_Human\_HCD\_20170816\_\*.msp" and a red-bordered "Edit filename pattern" button. At the bottom, there is a "Start downloading" button. Two yellow callout boxes provide instructions: the first points to the "Edit download URLs" button with the text "① 必要に応じて msp ファイルの取得先を編集" and the second points to the "Start downloading" button with the text "② 「Start downloading」ボタンを押し、ファイルの取得を開始".

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Filename pattern  
NIST\_Human\_HCD\_20170816\_\*.msp

Please,

1. copy the MSP file to C:/inetpub/mascot/sequence/NIST\_Human\_HCD\_20170816/current;
2. rename the file to match the filename pattern, e.g. NIST\_Human\_HCD\_20170816\_20170821.msp (in the pattern, \* matches any string of characters);
3. and refresh this page to update the database.

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
The screenshot shows the 'New spectral library definition' page in the Mascot Database Manager. The left sidebar has a 'Library' section with a 'Create new' button highlighted by a red box and labeled ①. The main form has a 'Library name (?)' field with 'NIST\_Human\_HCD\_20170816'. Below it, the 'Copy of (?)' radio button is selected and highlighted with a red box, with 'NIST\_Human\_HCD' selected in the dropdown menu. The 'Next' button is also highlighted with a red box. A yellow callout box labeled ② explains that 'Copy of' should be selected with 'NIST\_Human\_HCD' and the 'Next' button pushed.

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**Library name:**  
NIST\_Human\_HCD\_20170816

**Base directory (?)**  
C:/inetpub/mascot/sequence

Library files will be located in the subdirectory *NIST\_Human\_HCD\_20170816* of the base directory. The new directory will be created if it does not already exist.

**Existing NIST\_Human\_HCD files:**  
NIST\_Human\_HCD\_20160504.msp

**Copy files also (?)**  
 **Don't copy files (?)**

If you have chosen to copy files from the base definition, please do not close the browser or refresh this page after clicking Create. Copying may take some time if the files are large.

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The screenshot shows the Mascot Database Manager configuration page for the database 'NIST\_Human\_HCD\_20170816'. The page is divided into several sections:

- Database Manager:** Includes links for Databases (12), Parse rules (15), Scheduled updates (0), Running tasks (0), and Settings.
- Fasta:** Includes links for Enable predefined definition, Synchronise custom definitions, and Create new.
- Library:** Includes links for Enable predefined definition, Synchronise custom definitions, Create new, and Spectral library filters.

The main configuration area includes:

- Database: NIST\_Human\_HCD\_20170816** (with Copy and Delete buttons)
- Name:** NIST\_Human\_HCD\_20170816
- Database type:** Spectral library
- Database directory:** C:/inetpub/mascot/sequence/NIST\_Human\_HCD\_20170816/current (highlighted with a red box and labeled 'msp ファイル設置場所')
- Filename pattern:** NIST\_Human\_HCD\_20170816\_\*.msp (highlighted with a red box and labeled 'msp ファイルの名称パターン')
- No local MSP file:** Database files must be present before the database can be configured and activated. You have the following options.
  - Download from remote URL or copy from Mascot server hard disk (URLs are currently defined.)
  - Upload file using web browser
  - Copy file manually

A 'Next' button is located at the bottom of the configuration area.

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Filename pattern: NIST\_Human\_HCD\_20170816\_\*.msp

**Download files from URL or copy from Mascot Server disk**

FASTA file URL: ftp://chemdata.nist.gov/download/pep/.../human\_hcd\_selected.msp.tar.gz

**Edit download URLs**

Filename pattern: NIST\_Human\_HCD\_20170816\_\*.msp

**Edit filename pattern**

The downloaded files will be renamed...

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The screenshot shows the Mascot Database Manager interface. The browser address bar indicates the URL is localhost/mascot/x-cgi/db\_manager.pl?dbs.file-source=uplc. The main content area is titled 'Database: NIST\_Human\_HCD\_20170816' and includes fields for Name, Database type (Spectral library), Database directory, and Filename pattern. Below this is the 'Upload MSP file' section, which includes a note about large files and a 'Filename pattern' field. Two yellow callout boxes with red text provide instructions: ① 「参照」ボタンを押してアップロードするファイルを選び、 and ② 「Upload」ボタンを押します。 The '参照...' button and the 'Upload' button are circled in red in the screenshot.

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
The screenshot also shows the 'Database Status' page, which displays the status of various databases. The 'NIST\_Human\_HCD\_20170816' database is highlighted with a red box, and a yellow box highlights the text 'Database Status 画面で データベース構築確認'.

```
Name = PRIDE_Contaminants Family = C:/inetpub/mascot/sequence/PRIDE_Contaminants/current/PRIDE_Contaminants_*.msp
Filename = PRIDE_Contaminants_20160908.msp Pathname = C:/inetpub/mascot/sequence/PRIDE_Contaminants/current/PRIDE_Contaminants_20160908.msp
Status = In use
State Time = Wed Aug 16 09:57:56 # searches = 0
Mem mapped = NO Request to mem map = NO Request unmap = NO Mem locked = NO
Number of threads = 1 Current = YES Type = Spectral library

Name = SwissProt Family = C:/inetpub/mascot/sequence/SwissProt/current/SwissProt_*.fasta
Filename = SwissProt_2017_05.fasta Pathname = C:/inetpub/mascot/sequence/SwissProt/current/SwissProt_2017_05.fasta
Status = In use
State Time = Wed Aug 16 09:57:56 # searches = 0
Mem mapped = YES Request to mem map = YES Request unmap = YES Mem locked = YES
Number of threads = -1 Current = YES Type = Amino acid

Name = NIST_Human_HCD_20170816 Family = C:/inetpub/mascot/sequence/NIST_Human_HCD_20170816/current/NIST_Human_HCD_20170816_*.msp
Filename = NIST_Human_HCD_20170816_fixed.msp Pathname = C:/inetpub/mascot/sequence/NIST_Human_HCD_20170816/current/NIST_Human_HCD_20170816_fixed.msp
Status = Creating compressed files 12% complete
State Time = Mon Aug 21 11:29:11 # searches = 0
Mem mapped = NO Request to mem map = NO Request unmap = NO Mem locked = NO
Number of threads = 1 Current = NO Type = Spectral library
```

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
The screenshot shows the Mascot Database Manager interface. The top part displays the 'Database: OriginalSL' configuration page. The left sidebar contains navigation menus for 'Database Manager', 'Fasta', and 'Library'. The main content area shows the database name 'OriginalSL', its type 'Spectral library', and its directory 'C:/inetpub/mascot/sequence/OriginalSL/current'. The 'FASTA file source path' is set to 'C:/inetpub/mascot/sequence/PRIDE\_Contaminants/current/PRIDE\_Contaminants\_20160908.msp'. A red box highlights the 'Edit download URLs' button. A yellow callout box with a red arrow points to this button, containing the text: ① 「Edit download URLs」を押し、.

Below this, the 'Filename pattern' is 'OriginalSL\_\*.msp', with an 'Edit filename pattern' button. A note states: 'The downloaded files will be renamed to match the filename pattern.' A red box highlights the 'Start downloading' button. A yellow callout box with a red arrow points to this button, containing the text: ③ 元の画面に戻ったのちに「Start downloading」を押してダウンロード/コピーを実施する.

The bottom part of the screenshot shows the 'Library configuration: OriginalSL' page. The 'Main file URL or path to source file on Mascot Server hard disk' is the same path as above. A checkbox for 'Delete original file after updating' is present. A note states: 'The original file can only be deleted if Database Manager has sufficient permissions.' A red box highlights the 'Save' button. A yellow callout box with a red arrow points to this button, containing the text: ② 入力欄に、ダウンロードであれば URL を、コピーであれば path を記入して「Save」.

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The screenshot shows the Mascot Database Manager interface. The browser address bar indicates the URL: `logstos/mascot/x-cgi/db_manager.pl?dbs.file-source=upload&sub=dbs%3AOriginalSL`. The page title is "Database: OriginalSL".

On the left sidebar, there are sections for "Database Manager" (Databases (12), Parse rules (15), Scheduled updates (0), Running tasks (0), Settings), "Fasta" (Enable predefined definition, Synchronise custom definitions, Create new), and "Library" (Enable predefined definition, Synchronise custom definitions, Create new, Spectral library filters).

The main content area shows the following details for the "OriginalSL" database:

- Name:** OriginalSL
- Database type:** Spectral library
- Database directory:** C:/inetpub/mascot/sequence/OriginalSL/current
- Filename pattern:** OriginalSL\_\*.msp

The "Upload MSP file" section contains the following elements:

- Text: "Please select the MSP file to upload." and "Note that if the file is large (more than 50 megabytes), it is better to transfer the file to the Mascot Server hard disk by some other means than uploading via web browser."
- Filename pattern:** OriginalSL\_\*.msp (with an "Edit filename pattern" button)
- MSP file to upload:** A button labeled "ファイルを選択" (Select file) is highlighted with a red box. The text "選択されていません" (Not selected) is next to it.
- Text: "The uploaded file will be renamed to match the filename pattern."
- Buttons: "Previous" and "Upload" (highlighted with a red box).

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Name  
OriginalSL

Database type  
Spectral library

**Database directory**  
C:/inetpub/mascot/sequence/OriginalSL/current

Filename pattern  
OriginalSL\_\*.msp

### Copy files manually

Filename pattern  
OriginalSL\_\*.msp

Edit filename pattern

Please,

1. copy the MSP file to C:/inetpub/mascot/sequence/OriginalSL/current
2. rename the file to match the pattern, \* matches any string
3. and refresh this page to content

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

**Accession parse rule for OriginalSL**

Raw test data (10 samples)


```

1. sp|TRYP_PIG|
2. sp|GFP_AEQV1|
3. sp|K2M3_SHEEP|
4. sp|ALBU_BOVIN|
5. sp|K2C1_HUMAN|
6. sp|RS27A_HUMAN|
7. sp|MYG_HUMAN|
8. sp|EGF_HUMAN|
9. sp|FABPH_HUMAN|
10. sp|CATA_HUMAN|

```

Parse rule	Match	Extracted data
<input type="radio"/> ¥(.*)	10/10	1. sp TRYP_PIG  ▶9 more matches
<input type="radio"/> ¥([^\ ]*)	10/10	1. sp TRYP_PIG  ▶9 more matches

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## Accession parse rule for OriginalSL

### Raw test data (10 samples)

1. sp|TRYP\_PIG|
2. sp|GFP\_AEQVI|
3. sp|K2M3\_SHEEP|
4. sp|ALBU\_BOVIN|
5. sp|K2C1\_HUMAN|
6. sp|RS27A\_HUMAN|
7. sp|MYG\_HUMAN|
8. sp|EGF\_HUMAN|
9. sp|FABPH\_HUMAN|
10. sp|CATA\_HUMAN|

When editing a parse rule:

- Parse rules are Basic Regular Expressions, used to extract information from text strings. A full description can be found in *Appendix A: Basic Regular Expressions* in the *Mascot: Installation and Setup* manual.
- Most characters match themselves. The following characters are special, and have to be backslashed to match themselves: `.[¥*^$`
- The information to be extracted must be enclosed in backslashed parentheses. e.g. `¥(extract me¥)`
- A list of characters enclosed in square brackets means "match any from the list". If the first character in the list is a caret `^`, it means "match any except those on the list".
- An asterisk means "match as many as possible". Thus, `[^ ]*` means "match as many characters as possible up to a space".
- A period `.` represents any character.

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The screenshot shows the 'Description parse rule for OriginalSL' configuration page. It includes a list of raw test data, a table of suitable parse rules, and a section for unsuitable parse rules. The 'Choose' button is highlighted in red. Below the screenshot, a smaller window shows the 'MSP file' configuration screen, where the 'Next' button is also highlighted in red.

**Description parse rule for OriginalSL**

**Raw test data (10 samples)**

1. sp|TRYP\_PIG|
2. sp|GFP\_AEQV1|
3. sp|K2M3\_SHEEP|
4. sp|ALBU\_BOVIN|
5. sp|K2C1\_HUMAN|
6. sp|RS27A\_HUMAN|
7. sp|MYG\_HUMAN|
8. sp|ECGF\_HUMAN|
9. sp|FABPH\_HUMAN|
10. sp|CATA\_HUMAN|

**Suitable parse rules (3)**

Parse rule	Match	Extracted data
<input type="radio"/> .. ¥([^\] *¥)	10/10	1. TRYP_PIG ▶9 more matches
<input checked="" type="radio"/> ¥(. *¥)	10/10	1. sp TRYP_PIG  ▶9 more matches
<input type="radio"/> ¥([^\] *¥)	10/10	1. sp TRYP_PIG  ▶9 more matches

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

The screenshot shows the Mascot database manager interface. The main panel displays the configuration for the 'OriginalSL' database. The 'Database status' section shows 'Offline' and an 'Activate' button, which is highlighted with a red box. A yellow callout box points to this button with the text: 「Activate」 Database 構築開始.

Below the configuration page, a status log is visible. The log entry for 'OriginalSL' is highlighted with a red box. A yellow callout box points to the 'Status = In use' line with the text: 目的データベースの status 項目が「In use」となれば 使用可能.



```

Name = SwissProt Family = C:/inetpub/mascot/sequence/SwissProt/current
Filename = SwissProt_2017_05.fasta Pathname = C:/inetpub/mascot/sequence/SwissProt/current
Status = In use Status = In use
State Time = Tue Aug 22 16:50:56 # searches = 0
Mem mapped = YES Request to mem map = YES Request to mem map = YES
Number of threads = -1 Current = YES Type = Amino acid sequence

Name = OriginalSL Family = C:/inetpub/mascot/sequence/OriginalSL/current
Filename = OriginalSL_20170522.msp Pathname = C:/inetpub/mascot/sequence/OriginalSL/current
Status = Creating compressed files - 20% complete
State Time = Wed Aug 23 15:26:09 # searches = 0
Mem mapped = NO Request to mem map = NO Request to mem map = NO Mem locked = NO
Number of threads = 1 Current = NO Type = Spectral library

```

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The screenshot shows the 'Create spectral library from search results' page in the Mascot Database Manager. The form contains the following fields and options:

- Library name:** OriginalSL2
- Sequence directory:** C:/inetpub/mascot/sequence
- Reference database:** SwissProt
- Taxonomy:** Homo sapiens (human)
- Fragment mass tolerance:** 0.6 Da
- Peptide mass tolerance:** 500 ppm
- Buttons:** Previous, Create

Two yellow callout boxes provide instructions:

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The screenshot shows the Mascot Database Manager interface for the 'OriginalSL2' database. The left sidebar contains navigation menus for 'Database Manager', 'Fasta', and 'Library'. The main content area displays database details and a section titled 'Create MSP file from search results'. In this section, the 'Peptide match filters' are currently set to '(none)'. The 'Edit filters' button is highlighted with a red rectangular box. A yellow callout box with a red border and a circled '6' icon contains the Japanese text: 「Edit filters」ボタンを押す。

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The screenshot shows the 'Peptide match filters for OriginalSL2' configuration page. It provides instructions on how to set filters, such as 'The library must have at least one score or expect value filter, typically expect < 0.01.' Below the text, there is a form with an 'Expect value' input field, a dropdown menu set to '<math>\leq</math>', and a text box containing '0.01'. The 'Expect value' input field is highlighted with a red rectangular box. A yellow callout box with a red border and a circled '7' icon contains the Japanese text: Expect value など、取り込みに使用する条件を指定してください。

**Peptide match filters for OriginalSL2**

The library must have at least one score or expect value filter, typically `expect < 0.01`.

Each individual filter is in a filter group. To add more filters to the group, use the OR button. To add more groups, use the AND button. The peptide match must pass all filter groups to be accepted, but within each group, only one filter needs to succeed.

To remove a filter, leave its value field empty. To remove a filter group, remove all its filters.

Filters are used in two complementary ways:

1. When Database Manager chooses results files to process, only files that might contain suitable peptide matches are included.
2. When Database Manager loops over peptide matches in a results file, only matches that pass the filter are imported to the library.

For example, if you have a filter DB = SwissProt and no other DB filters, then only results files that were searched against SwissProt are processed. (Or in a multi-database search, had SwissProt as one of the databases.) When Database Manager loops over its peptide matches, only those that actually come from SwissProt are imported.

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The screenshot shows the Mascot Database Manager interface for the 'OriginalSL2' database. The sidebar on the left includes sections for 'Database Manager', 'Fasta', and 'Library'. The main content area displays database details such as Name, Database type, Database directory, and Filename pattern. A yellow callout box with a red border contains the text '㉑ 「Import search results」ボタンを押してください。' (㉑ Press the 'Import search results' button). The 'Import search results' button is highlighted with a red rectangle.

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The screenshot shows the 'Import search results in OriginalSL2' form. The 'Results file date range' field is highlighted with a red rectangle and contains the text 'From midnight (0:00) on 2017-07-03 to midnight (23:59) on 2017-07-06'. The 'Filepath wildcard' field contains the text '../data/\*/\*'. A yellow callout box with a red border contains the text '㉒ さらに、検索日による絞り込みを行い、' (㉒ Additionally, narrow down by search date). Another yellow callout box with a red border contains the text '㉓ 「Add import task to queue」ボタンを押してデータ取り込みを開始します。' (㉓ Press the 'Add import task to queue' button to start data import). The 'Add import task to queue' button is highlighted with a red rectangle.

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Database Manager

- Databases (12)
- Parse rules (17)
- Scheduled updates
- Running tasks (1)**
- Settings

**Running tasks**

This page will refresh every 10 seconds.

# Queued task

#	Queued task	Added
1	Scan for results files between 2017-07-03 and 2017-07-06 and import in OriginalSL2 dbman_scan_for_library.pl OriginalSL2 --start-date 2017-07-03 --end-date 2017-07-06 --filepath-wildcard ../data/*/*	Wed Aug 23 18:07:29 2017

Waiting for next queued task to start...

**Task log**

Process end time	Command
[Wed Aug 23 18:04:05 2017]	dbman_scan_for_library.pl OriginalSL2 --filepath-wildcard ../data/*/* --search --start-date 2017-08-07 --end-date 2017-08-1
[Mon Aug 14 16:23:17 2017]	dbman_download.pl NIST_Mouse_HCD
[Tue Jun 27 12:13:13 2017]	dbman_download.pl NIST_Human_HCD
[Wed May 31 15:10:18 2017]	dbman_scan_for_library.pl filters4 --search --start-date 2017-05-24 --end-date 2017-05-31

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Task processing is currently **enabled**. Task queue **control**

This page will refresh every 10 seconds.

No tasks in queue.

Title	Started	Progress
Scan for results files between 2017-07-03 and 2017-07-06 and import in OriginalSL2 dbman_scan_for_library.pl OriginalSL2 --start-date 2017-07-03 --end-date 2017-07-06 --filepath-wildcard ../data/*/*	Wed Aug 23 18:08:00 2017	<b>35.7</b> % ../data/F001281.dat: Argument "SwissProt" isn't numeric in sort at C:/inetpub/mascot/perl64/site/lib/DBManager/Workhorse/SLFilterCompiler.pm line 154.


**Task log**

Process end time	Command
[Wed Aug 23 18:04:05 2017]	dbman_scan_for_library.pl OriginalSL2 --filepath-wildcard ../data/*/* --search --start-date 2017-08-07 --end-date 2017-08-1

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
The screenshot shows the Mascot Database Manager interface. The main heading is "Import search results in MyLibrary01". Below this, there are several sections for configuration:

- Results file date range:** From midnight (0:00) on 1970-01-01 to
- Filepath wildcard:** The default is to look in the daily directory. The input field contains `../data/*/*`.
- Library:** Includes options for "Include files already imported" and "Delete existing library contents".

A yellow callout box contains the following Japanese text:

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The "Add import task to queue" button is highlighted with a red rectangle.

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