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- ☒ Web Server
 - 4 ☒ Common HTTP Features
 - ☒ Default Document
 - ☒ Directory Browsing
 - ☒ HTTP Errors
 - ☒ Static Content
 - ☐ HTTP Redirection
 - ☐ WebDAV Publishing
 - 4 ☒ Health and Diagnostics
 - ☒ HTTP Logging
 - ☐ Custom Logging
 - ☐ Logging Tools
 - ☐ ODBC Logging
 - ☐ Request Monitor
 - ☐ Tracing
 - 4 ☒ Performance
 - ☒ Static Content Compression
 - ☐ Dynamic Content Compression
 - 4 ☒ Security
 - ☒ Request Filtering
 - ☐ Basic Authentication
 - ☐ Centralized SSL Certificate Support
 - ☐ Client Certificate Mapping Authentication
 - ☐ Digest Authentication
 - ☐ IIS Client Certificate Mapping Authentication
 - ☐ IP and Domain Restrictions
 - ☐ URL Authorization
 - ☐ Windows Authentication
 - 4 ☒ Application Development
 - ☐ .NET Extensibility 3.5
 - ☐ .NET Extensibility 4.5
 - ☐ Application Initialization
 - ☐ ASP
 - ☐ ASP.NET 3.5
 - ☐ ASP.NET 4.5
 - ☒ CGI
 - ☐ ISAPI Extensions
 - ☐ ISAPI Filters
 - ☐ Server Side Includes
 - ☐ WebSocket Protocol
 - ☐ FTP Server
 - ☐ IIS Hostable Web Core
 - ☒ Management Tools
 - ☒ IIS Management Console
 - 4 ☒ IIS 6 Management Compatibility
 - ☒ IIS 6 Metabase Compatibility
 - ☒ IIS 6 Management Console
 - ☒ IIS 6 Scripting Tools
 - ☒ IIS 6 WMI Compatibility
 - ☒ IIS Management Scripts and Tools
 - ☐ Management Service

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Please type in your contact details to finish creating the licence.

- If any of the product and licence details shown below are incorrect, please contact Matrix Science support.
- Required fields are marked with an asterisk*.

The licence file will be sent to the email address you supply below. Make sure your spam filter does not reject this message by adding support@matrixscience.com to your [whitelist] or address book.

Product key

Product family Mascot Server

Licence type Permanent

Licenced features Version
Version
2 CPUs

Old product key

Email address*

Verify email address*

Full name*

Organisation*

Address line 1*

Address line 2

City*

Zip/postcode

State/region

Country*

Phone number

Fax number

If you are installing the software as an end user, please select one of the following options and your e-mail address here.

Email address (CC)

☒ Matrix Science may send me technical support information by email, such as the availability of new releases, upgrades, and support contracts.

☐ I do not wish to receive any communications from Matrix Science.

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Note: In either case, Matrix Science will not use your contact information to send unrelated advertising, and will not disclose your details to any third party.

MASCOT PROTEIN IDENTIFICATION SYSTEM

End-user Licence Agreement

IMPORTANT – PLEASE READ CAREFULLY: This End User Licence Agreement is a legally binding contract between you (either an individual or a single corporate entity) and Matrix Science Limited for the product identified above, which includes computer software, electronic documentation, any printed documentation, and any subsequent updates and supplements (the "Software").

By installing or using the Software, you agree to be bound by the terms of this agreement. If you do not agree to the terms of this agreement, we are unwilling to license the Software to you. In this case, do not install or use the Software. Return the package that included the Software to Matrix Science Limited or their authorised distributor within 30 days of receipt for a full refund.

1 Licence

Matrix Science Limited owns the copyright in the Software contained within this package and all other copies which you are authorised by this agreement to make.

This licence is personal to you (either an individual or a single corporate entity) as the purchaser of a

☐ I do not accept the terms of the licence agreement.



☒ I accept the terms of the licence agreement.*

Create Licence

Copyright © 1998-2020 Matrix Science Ltd. All Rights Reserved.

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Name = Human EST
Filename =
Status = Not in use
State Time = Thu Jan 1 09:00:00
Mem mapped = NO Request to mem
Number of threads = -1 Current

Name = NIST DSA IonTrap
Filename =
Status = Not in use
State Time = Thu Jan 1 09:00:00
Mem mapped = NO Request to mem
Number of threads = -1 Current

Mascot Server Licence Information

[Register a new product key](#) [View database status](#) [Reload this page](#)

Please include *all* the contents of this page when requested to provide this information to technical support.

Mascot Server version: 2.8.0.1
Licence path: C:\inetpub\mascot\config\licdb

Licence(s) found:

Product Key	Start	End	Status
Active:			
2021-08-05 *			OK
Feature: Mascot Server - Core functionality (v2.8)			
Feature: Mascot Server - CPU units (1)			
Feature: Mascot Server - Cores per CPU units (6)			
Company: mskk			
User: takaesu			
Distributor: Matrix Science KK			
Inactive:			
2021-08-05 *			OK
Feature: Mascot Server - Core functionality (v2.7)			
Feature: Mascot Server - CPU units (1)			
Feature: Mascot Server - Cores per CPU units (6)			
Company: mskk			
User: takaesu			
Distributor: Matrix Science KK			

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
Mascot Database Manager Configuration Import

Note:
Please disable other programs and scripts that may interfere with the Mascot Server.

Database Manager can keep your sequence data files.
If you choose to use Database Manager, you should keep your sequence data files.
If you do not wish to use Database Manager, you should delete your sequence data files.

Database Manager has two kinds of database definition:

1. Predefined database definition: Automatic (Synchronise). Synchronisation means the predefined (FASTA, reference, taxonomy) will be updated.
2. Custom database definition: Manually made.

You can choose to keep the existing definition (Synchronise). Synchronisation means the predefined (FASTA, reference, taxonomy) will be updated.

If you do not want to update the files now, choose Keep as Custom. You can synchronise at any time in the future by choosing Synchronise Custom Definitions from the Database menu.

Database definition	Action	Match quality
PRIDE_Contaminants	<input type="radio"/> Keep as Custom <input checked="" type="radio"/> Synchronise with PRIDE_Contaminants	perfect
SwissProt	<input type="radio"/> Keep as Custom <input checked="" type="radio"/> Synchronise with SwissProt_ID <input type="radio"/> Synchronise with SwissProt_AC	perfect poor Differences

Import database configuration by choosing Import.

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Latest predefined definitions files are from Thu Mar 8 20:27:25 2018 (FASTA databases: databases_1.xml) and Tue Nov 8 00:24:22 2016 (spectral libraries: libraries_1.xml).
Full database status is available on the [database status page](#).

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definition

Databases

Name	Status	New files	Tasks
IPI_human	* In use	Update	Deactivate 1 task in queue
NCBINr	* In use	Update	Deactivate 1 task in queue
OriginalDB	In use	(no URLs) Deactivate	
SwissProt	* In use	Update Deactivate	1 task in queue
cRAP	Offline	(no URLs) Activate	
EST_human	* Offline	Download No local FASTA file	1 task in queue

*) Entries marked with an asterisk are predefined definitions. Predefined definitions file was last updated Tue May 8 17:19:46 2012.

Full database status information is available on [the database status page](#).

[Refresh](#)

Proxy settings

You may need to define proxy settings to enable Database Manager to access the Internet, depending on your local network and Internet connection.

Proxy settings defined in:

- HTTP proxy strategy
- HTTPS proxy strategy
- FTP proxy strategy

HTTP test succeeds can access http://http-connectivity-test.matrixscience.com/db_manager.txt from the Mascot server machine.

HTTPS test succeeds can access https://https-connectivity-test.matrixscience.com/db_manager.txt from the Mascot server machine.

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Proxy settings

mascot.dat settings

The recommended proxy type is **Auto**, or failing that, **Environment**, **Registry**, **Specify**, or **WPAD**.

Proxy type

Specify

Host (HTTP proxy) (?)

Port (HTTP proxy) (?)

User (HTTP proxy) (?)

Password (HTTP proxy) (?)

☐ Use HTTP proxy for FTP (?)

☐ Use HTTP proxy for HTTPS (?)

Save

HTTP test **succeeds**; can access http://http-connectivity-test.matrixscience.com/db_ server machine.

HTTPS test **succeeds**; can access <https://https-connectivity-test.matrixscience.com/c> server machine.

(FTP proxy settings not tested.)

Back

Proxy type

- Auto
- Auto**
- None
- Environment
- Registry
- Specify
- WPAD

Auto means Database Manager will try to get proxy settings from the system registry.

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The screenshot shows the Mascot Database Manager interface. The left sidebar has two main sections: 'Fasta' and 'Library'. Under 'Fasta', the 'Enable predefined definition' link is highlighted with a yellow box. Under 'Library', the 'Enable predefined definition' link is also highlighted with a yellow box. The main content area displays a table of databases and spectral libraries. Two yellow callout boxes provide instructions: one for Fasta and one for Library.

Name	Mode	Type	Status	Latest task
PRIDE_Contaminants	predefined	SL	(waiting)	Get new files Deactivate
SwissProt	predefined	AA	In use	Get new files Deactivate

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Enable predefined definition

Library name: NIST_Human_IonTrap
Base directory (?): C:/inetpub/mascot/sequence
Library files will be located in the subdirectory NIST_Human_IonTrap directory. The new directory will be created if it does not exist.

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About to download or copy files

Database files need to be downloaded or copied before database configuration can be completed.
Download task is in queue position 1; waiting for downloading to start. See the task queue for more details. (Download task added Thu Feb 2 11:32:14 2017.)

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<対象生物種リスト:DB名XXXXに該当> Environmental , Fungi, Invertebrates, Mammals, Mus, Plants, Prokaryotes, Rodents, Vertebrates

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<対象生物種リスト:DB名YYYYに該当> arabidopsis, bovine, chicken , human, mouse, rat, zebrafish

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The screenshot displays the Mascot Database Manager interface. The top window shows the 'Database: NIST_Human_HCD' page with a sidebar containing 'Running tasks (1)'. The bottom window shows the 'Database Manager tasks' page, which includes a task list and a task log. A task 'Download files for NIST_Human_HCD' is shown in progress with a progress bar at 6.0%. The task log at the bottom shows a successful update for 'NIST_Human_IonTrap' with the message '(success) 'NIST_Human_IonTrap' successfully updated.'.

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(success) 'NIST_Human_IonTrap' successfully updated.

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